NKMAXBIO We support you, we believe in your research

Human VAMP-2 antibody

Catalog Number: AVA0710

PRODUCT INFORMATION

Catalog number

AVA0710

Clone No.

3E5

Product type

Monoclonal Antibody

UnitProt No.

P63027

NCBI Accession No.

NP 055047

Alternative Names

vesicle-associated membrane protein 2, VAMP2, SYB2, Synaptobrevin 2, VAMP-2, Synaptobrevin-2, Vesicle-associated membrane protein 2, vesicle-associated membrane protein 2 FLJ11460, RATVAMPB, RATVAMPIR, SYB, VAMP 2, Vesicle associated membrane protein 2, Vesicle-associated membrane protein 2 (synaptobrevin 2)

PRODUCT SPECIFICATION

Antibody Host

Mouse

Reacts With

Human, Mouse

Concentration

1mg/ml (determined by BCA assay)

Formulation

Liquid in. Phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% glycerol

Immunogen

Recombinant human Synaptobrevin 2 (1-89aa) purified from E. coli

Isotype

IgG1 kappa

Purification Note

By protein-G affinity chromatography

Application

ELISA, WB, ICC/IF

Usage

The antibody has been tested by ELISA, Western blot and ICC/IF analysis to assure specificity and reactivity. Since application varies, however, each investigation should be titrated by the reagent to obtain optimal results.



NKMAXBio We support you, we believe in your research

Human VAMP-2 antibody

Catalog Number: AVA0710

Storage

Can be stored at +2C to +8C for 1 week. For long term storage, aliquot and store at -20C to -80C. Avoid repeated freezing and thawing cycles.

BACKGROUND

Description


Synaptobrevin 2 (Vesicle-associated membrane, VAMP2), which is an 18 kDa integral membrane protein localized to the cytoplasmic surface of synaptic vesicle, consists of a proline-rich N-terminal region, a highly conserved hydrophilic domain, followed by a transmembrane anchor and a C-terminal tail. Synaptobrevin 2 is predominantly expressed in Langerhans islets and glomerular cells. The N-terminal domain of the protein forms a specific SNARE complex with the target membrane-associated t- or Q-SNAREs syntaxin 1 and SNAP-25. This antibody recognizes specifically synaptobrevin 2, but it also shows a low affinity to synaptobrevin 1.


General References

Andres C et al.,(1995) Proc Natl Acad Sci USA 92:5987-5991. Jacobsson G et al.,(1998) European Journal of Neuroscience 16:301-316. Scales SJ. et al (2002) J Biol. Chem. 227(31):28271-28279. Pennuto M et al., (2003) Mol Biol Cell. 14:4909-19.

DATA

Western blot analysis (WB)





The recombinant proteins (50ng) were resolved by SDS-PAGE, transferred to PVDF membrane and probed with anti-human VAMP-2 antibody (1:1000). Proteins were visualized using a goat anti-mouse secondary antibody conjugated to HRP and an ECL detection system.

Lane 1.: Recombinant human VAMP1 protein Lane 2.: Recombinant human VAMP2 protein Lane 3.: Recombinant human VAMP3 protein


The tissue lysate (50ug) was resolved by SDS-PAGE, transferred to PVDF membrane and probed with anti-human VAMP-2 antibody (1:1000). Proteins were visualized using a goat anti-mouse secondary antibody conjugated to HRP and an ECL detection system.

Lane 1.: Mouse brain Tissue lysate

NKMAXBio We support you, we believe in your research


Human VAMP-2 antibody

Catalog Number: AVA0710



The tissue lysate (20ug) was resolved by SDS-PAGE, transferred to PVDF membrane and probed with anti-human VAMP-2 antibody (1:2000). Proteins were visualized using a goat anti-mouse secondary antibody conjugated to HRP and an ECL detection system. Lane 1.: Rat brain Tissue lysate

Immunocytochemistry/Immunofluorescence (ICC/IF)



ICC/IF analysis of VAMP-2 in U87MG cells. The cell was stained with AVA0710 (1:100). The secondary antibody (green) was used Alexa Fluor 488. DAPI was stained the cell nucleus (blue).

