NKMAXBIO We support you, we believe in your research

Recombinant rat TNF-alpha protein

Catalog Number: ATGP3128

PRODUCT INFORMATION

Expression system

E.coli

Domain

80-235aa

UniProt No.

P16599

NCBI Accession No.

NP 036807

Alternative Names

Tumor necrosis factor, Cachectin, TNF-alpha, Tumor necrosis factor ligand superfamily member 2, TNFSF2, TNF-a, TNFA, TNF

PRODUCT SPECIFICATION

Molecular Weight

19.9 kDa (181aa) confirmed by MALDI-TOF

Concentration

1mg/ml (determined by absorbance at 280nm)

Formulation

Liquid in. Phosphate-Buffered Saline (pH 7.4) containing 10% glycerol, 1mM DTT

Purity

> 90% by SDS-PAGE

Tag

His-Tag

Application

SDS-PAGE

Storage Condition

Can be stored at +2C to +8C for 1 week. For long term storage, aliquot and store at -20C to -80C. Avoid repeated freezing and thawing cycles.

BACKGROUND

Description

The also known as Tumor necrosis factor. The Protein acts as a cytokine binds TNF receptors plays a role in regulation of cell proliferation, induction of apoptosis, and inflammatory response. Recombinant rat Tnf, fused to His-tag at N-terminus, was expressed in E. coli and purified by using conventional chromatography techniques.

Amino acid Sequence

MGSSHHHHHH SSGLVPRGSH MGSHMLRSSS ONSSDKPVAH VVANHOAEEO LEWLSORANA LLANGMDLKD



NKMAXBio We support you, we believe in your research

Recombinant rat TNF-alpha protein

Catalog Number: ATGP3128


NQLVVPADGL YLIYSQVLFK GQGCPDYVLL THTVSRFAIS YQEKVSLLSA IKSPCPKDTP EGAELKPWYE PMYLGGVFQL EKGDLLSAEV NLPKYLDITE SGQVYFGVIA L

General References

Furuya T., et al. (2001) Genes Immun. 2:229-232. Warle M.C., et al. (2005) Transpl. Immunol. 14:77-82.

DATA

SDS-PAGE



15% SDS-PAGE (3ug)

3ug by SDS-PAGE under reducing condition and visualized by coomassie blue stain.

