NKMAXBIO We support you, we believe in your research

Murine NK Cell Activity ELISA Test Kit

Catalog Number: ATGBF001

PRODUCT INFORMATION

Catalog number

ATGBF001

PRODUCT SPECIFICATION

Description

The Murine NK Activity Kit is an IFN-y quantitation assay for plasma samples collected and prepared with the NK Cell activating agent. Murine NK Activity Kit is intended for research use only, for the measuring of Natural Killer Cell Activity (NKA) in Mice. **Kit insert English version download link**

Storage

Activator (NK Cell activating agent): -20°C / Other components: 2-8°C

Product type

ELISA assay

