NKMAXBio We support you, we believe in your research

MERS-CoV Spike RBD antibody

Catalog Number: ATGA0602

PRODUCT INFORMATION

Catalog number

ATGA0602

Clone No.

AT18D10

Product type

Monoclonal antibody

UnitProt No.

K0BRG7

NCBI Accession No.

AFS88936

Alternative Names

Middle East respiratory syndrome coronavirus, Human betacoronavirus 2c EMC/2012, MERS-CoV, MERS, MERS-CoV RBD, MERS RBD, receptor binding domain, RBD, Spike RBD protein

PRODUCT SPECIFICATION

Antibody Host

Mouse

Reacts With

MERS-CoV

Concentration

1mg/ml (determined by BCA assay)

Formulation

Liquid in. Phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% glycerol

Immunogen

Recombinant MERS-CoV Spike RBD (358-606aa) purified from Baculovirus

Isotype

IgG1 kappa

Purification Note

By protein-A affinity chromatography

Application

ELISA, WB

Usage

The antibody has been tested by ELISA analysis to assure specificity and reactivity. Since application varies, however, each investigation should be titrated by the reagent to obtain optimal results.



NKMAXBIO We support you, we believe in your research

MERS-CoV Spike RBD antibody

Catalog Number: ATGA0602

Storage

Can be stored at +2C to +8C for 1 week. For long term storage, aliquot and store at -20C to -80C. Avoid repeated freezing and thawing cycles.

BACKGROUND

Description


MERS-CoV, which causes the Middles East Respiratory Syndrome (MERS), belongs to a family of viruses known as coronaviruses. MERS-CoV was first identified in the Kingdom of Saudi Arabia in 2012, which is a single and positive stranded RNA virus. Dromedary camels are widely considered as the source of the transmission of MERS-CoV. The rate of human transmission among household contacts of MERS patients has been approximately 5 % based on serological analysis. MERS-CoV has four structural proteins, known as the S (spike), E (envelope), M (membrane), and N (nucleocapsid) proteins. The spike protein, responsible for allowing the virus to attach to and fuse with the membrane of a host cell and is a large type I transmembrane protein containing two subunits, S1 and S2. S1 mainly contains a receptor binding domain (RBD), which is responsible for recognizing the cell surface receptor. S2 contains basic elements needed for the membrane fusion. MERS-CoV S mediates viral attachment and fusion to human cells via human cellular receptor DPP4, also known as CD26. The S protein plays key parts in the induction of neutralizing-antibody and T-cell responses, as well as protective immunity.

General References

Junghyun Goo., et al. (2020) Virus Res. 278:197863. Yan-Hua Li., et al. (2019) Engineering. 5:940-947. Lingshu Wang., et al. (2018) J Virol. 92:e02002-2017. Nicolas Papageorgiou., et al. (2016) Acta Crystallogr D Struct Biol. 72:192-202. Xiao-Yan Che., et al. (2004) J Clin Microbiol. 42:2629-2635.

DATA

Western blot analysis (WB)



Additional Information

The recombinant proteins (50ng) were resolved by SDS-PAGE, transferred to PVDF membrane and probed with anti-MERS Spike RBD antibody (1:1000). Proteins were visualized using a goat anti-mouse secondary antibody conjugated to HRP and an ECL detection system.

Lane 1.: Recombinant MERS Spike protein (ATGP3979)

Lane 2.: Recombinant MERS Spike S1 protein (ATGP3980)

Lane 3.: Recombinant MERS Spike S2 protein (ATGP3981)


Lane 4.: Recombinant MERS Spike RBD protein (ATGP3982)



NKMAXBio We support you, we believe in your research

MERS-CoV Spike RBD antibody

Catalog Number: ATGA0602



ELISA: MERS Spike RBD Antibody (1ug/ml) specifically recognizes MERS Spike RBD recombinant protein, but not interacted MERS Spike, Spike S1 and Spike S2 recombinant protein recombinant protein in ELISA.

