NKMAXBIO We support you, we believe in your research

Human PSMD11 antibody

Catalog Number: ATGA0446

PRODUCT INFORMATION

Catalog number

ATGA0446

Clone No.

AT2C7

Product type

Monoclonal Antibody

UnitProt No.

000231

NCBI Accession No.

NP 002806

Alternative Names

proteasome 26S non-ATPase subunit 11, S9, Rpn6, p44.5, MGC3844

PRODUCT SPECIFICATION

Antibody Host

Mouse

Reacts With

Human

Concentration

1mg/ml (determined by BCA assay)

Formulation

Liquid in. Phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% glycerol

Immunogen

Recombinant human PSMD11 (1-422aa) purified from E. coli

Isotype

IgG2b kappa

Purification Note

By protein-A affinity chromatography

Application

ELISA, WB, ICC/IF

Usage

The antibody has been tested by ELISA, Western blot analysis and ICC/IF to assure specificity and reactivity. Since application varies, however, each investigation should be titrated by the reagent to obtain optimal results.

Storage



NKMAXBIO We support you, we believe in your research

Human PSMD11 antibody

Catalog Number: ATGA0446

Can be stored at +2C to +8C for 1 week. For long term storage, aliquot and store at -20C to -80C. Avoid repeated freezing and thawing cycles.

BACKGROUND

Description


The 26S proteasome is a multicatalytic proteinase complex with a highly ordered structure composed of 2 complexes, a 20S core and a 19S regulator. The 20S core is composed of 4 rings of 28 non-identical subunits; 2 rings are composed of 7 alpha subunits and 2 rings are composed of 7 beta subunits. The 19S regulator is composed of a base, which contains 6 ATPase subunits and 2 non-ATPase subunits, and a lid, which contains up to 10 non-ATPase subunits. PSMD11 is a member of the proteasome subunit 59 family that functions as a nonATPase subunit of the 19S regulator and is phosphorylated by AMP-activated protein kinase.

General References

Ogris E., et al. (1999) J. Biol. Chem. 274: 14382-14391. Gagnon S.N., et al. (2002) Biochem. J. 368: 263-271.

DATA


Western blot analysis (WB)



The cell lysates (40ug) were resolved by SDS-PAGE, transferred to PVDF membrane and probed with anti-human PSMD11 antibody (1:1000). Proteins were visualized using a goat anti-mouse secondary antibody conjugated to HRP and an ECL detection system.

Lane 1.: HeLa cell lysate Lane 2.: 293T cell lysate Lane 3.: U87-MG cell lysate Lane 4.: NIH-3T3 cell lysate Lane 5.: PC3 cell lysate Lane 6.: TF1 cell lysate Lane 7.: MCF7 cell lysate

Immunocytochemistry/Immunofluorescence (ICC/IF)



ICC/IF analysis of PSMD11 in HeLa cells. The cell was stained with ATGA0446 (1:100). The secondary antibody (green) was used Alexa Fluor 488. DAPI was stained the cell nucleus (blue).

