

# Human CNDP2/CPGL antibody

Catalog Number: ATGA0562

## PRODUCT INFORMATION

---

**Catalog number**

ATGA0562

**Clone No.**

AT15E5

**Product type**

Monoclonal antibody

**UnitProt No.**

Q96KP4

**NCBI Accession No.**

NP\_060705

**Alternative Names**

Cytosolic non-specific dipeptidase, CN2, CPGL, HsT2298, PEPA

**Additional Information**

This product was produced from tissue culture supernatant.

## PRODUCT SPECIFICATION

---

**Antibody Host**

Mouse

**Reacts With**

Human

**Concentration**

1mg/ml (determined by BCA assay)

**Formulation**

Liquid in. Phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% glycerol

**Immunogen**

Recombinant human CNDP2 (1-475aa) purified from E. coli

**Isotype**

IgG1 kappa

**Purification Note**

By protein-A affinity chromatography

**Application**

ELISA, WB, ICC/IF

**Usage**

The antibody has been tested by ELISA, Western blot and ICC/IF analysis to assure specificity and reactivity. Since application varies, however, each investigation should be titrated by the reagent to obtain optimal results.

# Human CNDP2/CPGL antibody

Catalog Number: ATGA0562

## Storage

Can be stored at +2C to +8C for 1 week. For long term storage, aliquot and store at -20C to -80C. Avoid repeated freezing and thawing cycles.

## BACKGROUND

### Description


CNDP2, known as tissue carnosinase and peptidase A, is a nonspecific dipeptidase rather than a selective carnosinase and is a secreted peptidase homologous to M20 peptidases. Also known as glutamate carboxypeptidase-like protein 1 (CPGL), CNDP2 is ubiquitously expressed throughout human tissue. CNDP2 Hydrolyzes a variety of dipeptides including L-carnosine but has a strong preference for Cys-Gly and may be play a role as tumor suppressor in hepatocellular carcinoma (HCC) cells.

### General References

Zhang, P., et al. (2006) Cancer res 12: 6617-6625.  
Kaur, H., et al. (2009) J.Biol.Chem 284: 14493-14502.  
Licker, V., et al. (2012) J proteomics 375(15): 4656-67.

## DATA

### Western blot analysis (WB)




The tissue lysate (40ug) was resolved by SDS-PAGE, transferred to PVDF membrane and probed with anti-human CNDP2/CPGL antibody (1:1000). Proteins were visualized using a goat anti-mouse secondary antibody conjugated to HRP and an ECL detection system.  
Lane 1.: Mouse kidney tissue lysate

### Immunocytochemistry/Immunofluorescence (ICC/IF)

## Human CNDP2/CPGL antibody

Catalog Number: ATGA0562



ICC/IF analysis of CNDP2/CPGL in HeLa cells. The cell was stained with ATGA0562 (1:100). The secondary antibody (green) was used Alexa Fluor 488. DAPI was stained the cell nucleus (blue).