

SNRPF cDNA

Catalog Number: ATGD0105

PRODUCT INFORMATION

Catalog number

ATGD0105

Product type

cDNA

Species

Human

NCBI Accession No.

NP_003086.1

Alternative Names

Sm-F, SMF, snRNP-F

mRNA Refseq

NM_003095.2

OMIM

603541

Chromosome location

12q23.1

PRODUCT SPECIFICATION

Formulation

Lyophilized

Storage

Store the plasmid at -20C.

cDNA Size

261bp

Preparation before usage

1. Centrifuge at 7000rpm for 1 minute.
 2. Carefully open the vial and add 100ul of sterile water to dissolve the DNA.
- Each tube contains approximately 10ug of lyophilized plasmid.

Vector description

This shuttle vector contains the complete ORF. It is inseted BamH I to Xho I. The gene insert contains multiple cloning sites which can be used to easily cut and transfer the gene and recombination site into your expression vector.

Cloning Vector

pATGen (puc19-derived cloning vector)

General Description

SNRPF cDNA

Catalog Number: ATGD0105

DATA

Sequence nucleotides

ATGAGTTTACCCCTCAATCCCAAACCTTTCCTCAATGGACTAACAGGAAAGCCAGTGATGGTGAAACTTAAGTGGGGAATGG
AGTACAAGGGCTATCTGGTATCTGTAGATGGCTACATGAACATGCAGCTTGCAAATACAGAAGAATACATAGATGGAGCTTT
GTCTGGACATCTGGGTGAAGTTTTAATAAGGTGTAATAATGTCCTTTATATCAGAGGTGTGGAAGAAGAGGAAGAAGATGG
GGAAATGAGAGAATAG

Transaction Sequence

MSLPLNPKPF LNGLTGKPM VKLKWGMEYK GYLVSDGYM NMQLANTEY IDGALSGHLG EVLIRCNNVL YIRGVEEEE
DGEMRE