

ST20 cDNA

Catalog Number: ATGD0269

PRODUCT INFORMATION

Catalog number

ATGD0269

Product type

cDNA

Species

Human

NCBI Accession No.

NP_001094349.1

Alternative Names

HCCS-1

mRNA Refseq

NM_001100879.1

OMIM

Chromosome location

15q25.1

PRODUCT SPECIFICATION

Formulation

Lyophilized

Storage

Store the plasmid at -20C.

cDNA Size

240bp

Preparation before usage

1. Centrifuge at 7000rpm for 1 minute.
2. Carefully open the vial and add 100ul of sterile water to dissolve the DNA. Each tube contains approximately 10ug of lyophilized plasmid.

Vector description

Cloning Vector

pATGen (puc19-derived cloning vector)

General Description

DATA

Sequence nucleotides

ST20 cDNA

Catalog Number: ATGD0269

ATGGCGGATCTCGGCTCACTGCAACCTCTGTCTCCAGGTTTCAGGAAAATGGCTTTGTAAAGAAGCTTGAGCCTAAATCTG
GCTGGATGACTTTTCTAGAAGTTACAGGAAAGATCTGTGAAATGCTCTTCTGTCCTGAAGCAATACTGTTGACCAGAAAGGA
CACTCCATATTGTGAAACCGGCCTAATTTTTCTGACTCTTACGAAAACGATTGCCAACACATACTTCTACTTTTAA

Transaction Sequence

MARSRLTATS VSQVQENGFV KKLEPKSGWM TFLEVTGKIC EMLFCPEAIL LTRKDTPYCE TGLIFLTLTK TIANITYFYF