

TTC32 cDNA

Catalog Number: ATGD0272

PRODUCT INFORMATION

Catalog number

ATGD0272

Product type

cDNA

Species

Human

NCBI Accession No.

NP_001008238.1

Alternative Names

mRNA Refseq

NM_001008237.1

OMIM

Chromosome location

2p24.1

PRODUCT SPECIFICATION

Formulation

Lyophilized

Storage

Store the plasmid at -20C.

cDNA Size

456bp

Preparation before usage

1. Centrifuge at 7000rpm for 1 minute.
2. Carefully open the vial and add 100ul of sterile water to dissolve the DNA. Each tube contains approximately 10ug of lyophilized plasmid.

Vector description

This shuttle vector contains the complete ORF. It is inseted BamH I to Xho I. The gene insert contains multiple cloning sites which can be used to easily cut and transfer the gene and recombination site into your expression vector.

Cloning Vector

pATGen (puc19-derived cloning vector)

General Description

TTC32 cDNA

Catalog Number: ATGD0272

DATA

Sequence nucleotides

ATGGAAGGACAGCGGCAAGAAAGCCACGCAACCCTAACACTCGCCCAGGCTCATTCAACAATGGAGAGTACGCGGAGGC
CGAGGCACTGTACTCCGCTTACATTCGCCGGTGCGCTTGC GCGGCCTCCAGCGACGAGAGTCCCGGGAGCAAATGCAGCC
CTGAGGATTTGGCTACTGCATATAACAACAGGGGGCAAATCAAGTACTTCAGGGTTGATTTTTATGAAGCCATGGATGACTA
CACATCTGCCATAGAAGTCCAACCCAATTTTGAAGTTCCATATTACAACAGAGGGTTGATACTGTATAGGCTGGGATATTTTG
ATGATGCTTTGGAAGATTTCAAGAAGGTCTTAGACTTAAATCCTGGATTTCAAGATGCTACTTTGAGCTTAAAACAGACTATT
CTAGACAAAGAAGAAAAACAAAGAAGAAATGTTGCAAAAAATTATTGA

Transaction Sequence

MEGQRQESHA TLTLAQAHFN NGEYAEAEAL YSAYIRRCAC AASSDESPGS KCSPEDLATA YNNRGQIKYF RVDFYEAMDD
Y TSAIEVQPN FEVPYYNRGL ILYRLGYFDD ALEDFKKVLD LNPGFQDATL SLKQTILDKE EKQRRNVAKN Y