

## EID2B cDNA

Catalog Number: ATGD0282

### PRODUCT INFORMATION

---

**Catalog number**

ATGD0282

**Product type**

cDNA

**Species**

Human

**NCBI Accession No.**

NP\_689574.1

**Alternative Names**

EID-3

**mRNA Refseq**

NM\_152361.1

**OMIM**

**Chromosome location**

19q13.2

### PRODUCT SPECIFICATION

---

**Formulation**

Lyophilized

**Storage**

Store the plasmid at -20C.

**cDNA Size**

486bp

**Preparation before usage**

1. Centrifuge at 7000rpm for 1 minute.
2. Carefully open the vial and add 100ul of sterile water to dissolve the DNA. Each tube contains approximately 10ug of lyophilized plasmid.

**Vector description**

This shuttle vector contains the complete ORF. It is inseted Nde I to Xho I. The gene insert contains multiple cloning sites which can be used to easily cut and transfer the gene and recombination site into your expression vector.

**Cloning Vector**

pATGen (puc19-derived cloning vector)

**General Description**

## EID2B cDNA

Catalog Number: ATGD0282

### DATA

---

#### Sequence nucleotides

ATGGCGGAGCCGACTGGGCTGTTGGAGATGTCCGAGCTCCCCGGAGATAGCAGTGTCCCACAGGTGGGCACCGCGAGTG  
GCGTCAGCGACGTACTGCGGGGGGCAGTCGGCGGGCGGGGTTTCGGGTGCAGGAGGCCCGGGAAGGCCAGTGGCCGAAG  
CTGCGCGGTCCATGGCGCGGATGCCGGGCCCTGTGCCCGGGCCCATCCCCAGCAGCGTCCCGGGCCTGGCCTCCGCGCC  
AGACCCCATCAGCAGCTCGCTTTCTTAGAAATTAACCGGCAGCTGTTGTTCCGCGAGTATCTGGATGGTAGCTCCATGATT  
CCGGTCAGATTACTACGGGATTTTCGAGGAACGCCGCGAGGCTGTTTGTGGAAGGCTGCAAGGCGAGGGAAGCAGCCTTTGA  
CGCGGATCCCCCGCAGATGGACTTCGCTGCTGTCGCGTTTACGGTAGCGCTGACTGCCTCCGAGGCCCTCAGTCCTCTGGC  
CGACTGA

#### Transaction Sequence

MAEPTGLLEM SELPGDSSVP QVGTASGVSD VLRGAVGGGV RVQEAREGPV AEAARSMARM PGPVPGPIPS SVPGLASAPD  
PHQQLAFLEI NRQLLFREYL DGSSMIPVRL LRDFEERRRL FVEGCKAREA AFDADPPQMD FAAVAFTVAL TASEALSPLA D