NKMAXBio We support you, we believe in your research

SF3B5 cDNA

Catalog Number: ATGD0304

PRODUCT INFORMATION

Catalog number

ATGD0304

Product type

cDNA

Species

Human

NCBI Accession No.

NP 112577.1

Alternative Names

SF3b10, Ysf3

mRNA Refseq

NM_031287.2

OMIM

Chromosome location

6q24.2

PRODUCT SPECIFICATION

Formulation

Lyophilized

Storage

Store the plasmid at -20C.

cDNA Size

261bp

Preparation before usage

- 1. Centrifuge at 7000rpm for 1 minute.
- $2. \ \mbox{Carefully open the vial and add 100ul of sterile water to dissolve the DNA.}$

Each tube contains approximately 10ug of lyophilized plasmid.

Vector description

This shuttle vector contains the complete ORF. It is inseted BamH I to Xho I. The gene insert contains multiple cloning sites which can be used to easily cut and transfer the gene and recombination site into your expression vector.

Cloning Vector

pATGen (puc19-derived cloning vector)

General Description



NKMAXBio We support you, we believe in your research

SF3B5 cDNA

Catalog Number: ATGD0304

DATA

Sequence nucleotides

ATGACTGACCGCTACACCATCCATAGCCAGCTGGAGCACCTGCAGTCCAAGTACATCGGCACGGGCCACGCCGACACCACC AAGTGGGAGTGGCTGGAACCAACACCGCGACTCGTACTGCTCCTACATGGGCCACTTCGACCTTCTCAACTACTTCGCC ATTGCGGAGAATGAGAGCAAAGCGCGAGTCCGCTTCAACTTGATGGAAAAGATGCTTCAGCCTTGTGGACCGCCAGCCGAC AAGCCCGAGGAGAACTGA

Transaction Sequence

MTDRYTIHSQ LEHLQSKYIG TGHADTTKWE WLVNQHRDSY CSYMGHFDLL NYFAIAENES KARVRFNLME KMLQPCGPPA DKPEEN

