

NKMAX Protein Marker (7 bands, 8.5-70 kDa)

Catalog Number: PMR0401

NKMAX Protein Marker (7 bands, 8.5-70 kDa)

Catalog Number: PMR0401

PRODUCT INFORMATION

Expression system

E.coli

UniProt No.

N/A

Additional Information

Use 5ul per lane and you need not to heat before use

PRODUCT SPECIFICATION

Formulation

Liquid in. 12mM Tris-HCl buffer (pH 6.8) 5% glycerol, 0.4% SDS 2.88mM 2-mercaptoethanol, 0.02% bromophenol blue

Tag

Non-Tagged

Application

SDS-PAGE

Storage Condition

Can be stored at +2C to +8C for 1 week. For long term storage, aliquot and store at -20C to -80C. Avoid repeated freezing and thawing cycles.


BACKGROUND

Description

ATGen Mid-Range Protein Marker-1 consists of Seven precisely sized recombinant proteins (M. W. 70, 57, 40, 28, 18, 13.5, 8.5 kDa). The concentration of each protein is 0.5 ug/ul.

DATA

SDS-PAGE



15% SDS-PAGE (5ul)

3ug by SDS-PAGE under reducing condition and visualized by coomassie blue stain.