

Recombinant E.coli Carbonic anhydrase 2/CA2 protein

Catalog Number: CAN0801

PRODUCT INFORMATION

Expression system

E.coli

Domain

1-220aa

UniProt No.

P61517

NCBI Accession No.

NP_414668

Alternative Names

Carbonate dehydratase, CAN, ECK0125, JW0122, yadF

PRODUCT SPECIFICATION

Molecular Weight

27.2 kDa (240aa) confirmed by MALDI-TOF

Concentration

1mg/ml (determined by Bradford assay)

Formulation

Liquid in. 20mM Tris-HCl buffer (pH 8.0) containing 1mM DTT, 10% glycerol

Purity

> 95% by SDS-PAGE

Tag

His-Tag

Application

SDS-PAGE

Storage Condition

Can be stored at +2C to +8C for 1 week. For long term storage, aliquot and store at -20C to -80C. Avoid repeated freezing and thawing cycles.

BACKGROUND

Description

Carbonic anhydrase (CA) is an enzyme that catalyses rapid conversion of carbon dioxide to bicarbonate and protons ($\text{CO}_2 + \text{H}_2\text{O} \rightleftharpoons \text{HCO}_3^- + \text{H}^+$). Most carbonic anhydrases contain a zinc ion in their active site and the primary function of this enzyme is known to maintain acid-base balance in blood and other tissues, and to help transport carbon dioxide of tissues. Carbonic anhydrases have been found in all kingdoms of life. Recombinant carbonic anhydrase fused to His-tag, was expressed in E. coli and purified by conventional chromatography techniques.

Recombinant E.coli Carbonic anhydrase 2/CA2 protein

Catalog Number: CAN0801

Amino acid Sequence


<MGSSHHHHHH SSGLVPRGSH> MKDIDTLISN NALWSKMLVE EDPGFFEKLA QAQKPRFLWI GCSDSRVPAE
RLTGLEPGEL FVHRNVANLV IHTDLNCLSV VQYAVDVLEV EHHICGHYG CGGVQAAVEN PELGLINNWL LHIRDIWFKH
SSLLGEMPQE RRLDTLCELN VMEQVYNLGH STIMQSAWKR GQKVTIHGWA YGIHDGLLRD LDVTATNRET LEQRYRHGIS
NLKLLKHANHK

General References

Lindskog S., et al :(1997) Pharmacol Ther.74(1):1-20.
Sawaya MR., et al (2006) J Biol Chem. 281(11):7546-55

DATA

SDS-PAGE



3ug by SDS-PAGE under reducing condition and visualized by coomassie blue stain