NKMAXBIO We support you, we believe in your research

Recombinant human S100A4 protein

Catalog Number: ATGP0319

PRODUCT INFORMATION

Expression system

E.coli

Domain

1-101aa

UniProt No.

P26447

NCBI Accession No.

NP 002952

Alternative Names

S100 calcium-binding protein A4, S100 calcium-binding protein A4, 18A2, 42A, CAPL, FSP1, MTS1, P9KA, PEL98, Metastasin, S100A4, S100 calcium-binding protein A4 Protein S100 A4, S100 calcium binding protein A4.

PRODUCT SPECIFICATION

Molecular Weight

13.8 kDa (121aa) confirmed by MALDI-TOF

Concentration

1mg/ml (determined by Bradford assay)

Formulation

Liquid in. 20mM Tris-HCl buffer (pH 8.0) containing 10% glycerol

Purity

> 95% by SDS-PAGE

Tag

His-Tag

Application

SDS-PAGE

Storage Condition

Can be stored at +2C to +8C for 1 week. For long term storage, aliquot and store at -20C to -80C. Avoid repeated freezing and thawing cycles.

BACKGROUND

Description

S100A4, also known as S100 calcium-binding protein A4, belongs to the S100 super-family of proteins containing 2 EF-hand calcium binding domains. This protein is ubiquitously overexpressed and is localized in the cytoplasm and/or nucleus. S100A4 may play a role in motility, invasion, and tubulin polymerization. Chromosomal rearrangements and altered expression of the S100A4 gene have been implicated in tumor metastasis. Recombinant human S100A4 protein, fused to His-tag at N-terminus, was expressed in E. coli and purified by



NKMAXBio We support you, we believe in your research

Recombinant human S100A4 protein

Catalog Number: ATGP0319

using conventional chromatography.

Amino acid Sequence


MGSSHHHHHH SSGLVPRGSH MACPLEKALD VMVSTFHKYS GKEGDKFKLN KSELKELLTR ELPSFLGKRT DEAAFQKLMS NLDSNRDNEV DFQEYCVFLS CIAMMCNEFF EGFPDKQPRK K

General References

Sherbet GV., et al. (2009) Cancer Lett. 280(1):15-30. Ambartsumian N., et al. (1995) Gene. 159(1):125-30.

DATA

SDS-PAGE



3ug by SDS-PAGE under reducing condition and visualized by coomassie blue stain.

