NKMAXBIO We support you, we believe in your research

Recombinant human ERK2/MAPK1 protein

Catalog Number: ATGP0445

PRODUCT INFORMATION

Expression system

E.coli

Domain

1-360aa

UniProt No.

P28482

NCBI Accession No.

NP 620407

Alternative Names

PRKM2, PRKM1, p42MAPK, p41MAPK, p41, p40, p38, Mitogen-activated protein kinase 1 MAP kinase 2, Mitogen-activated protein kinase 1, MAPK2, MAPK1, Extracellular Signal Regulated Kinase 2, ERT1, ERK2, ERK, ,

PRODUCT SPECIFICATION

Molecular Weight

43.5 kDa (380aa) confirmed by MALDI-TOF

Concentration

1mg/ml (determined by Bradford assay)

Formulation

Liquid in. 20mM Tris-HCl buffer (pH 8.0) containing 10% glycerol, 1mM DTT, 0.1 M NaCl.

Purity

> 95% by SDS-PAGE

Tag

His-Tag

Application

SDS-PAGE

Storage Condition

Can be stored at +2C to +8C for 1 week. For long term storage, aliquot and store at -20C to -80C. Avoid repeated freezing and thawing cycles.

BACKGROUND

Description

MAPK1, also known as ERK (extracellular signal-regulated kinase), acts as an integration point for multiple biochemical signals, and is involved in a wide variety of cellular processes such as proliferation, differentiation, transcription regulation and development. The activation of this kinase requires its phosphorylation by upstream kinases. upon activation, this kinase translocates to the nucleus of the stimulated cells, where it phosphorylates nuclear targets. Recombinant human MAPK1 protein, fused to His-tag at N-terminus, was expressed in E. coli and



NKMAXBio We support you, we believe in your research

Recombinant human ERK2/MAPK1 protein

Catalog Number: ATGP0445

purified by using conventional chromatography techniques.

Amino acid Sequence


MGSSHHHHHH SSGLVPRGSH MAAAAAAGAG PEMVRGQVFD VGPRYTNLSY IGEGAYGMVC SAYDNVNKVR VAIKKISPFE HQTYCQRTLR EIKILLRFRH ENIIGINDII RAPTIEQMKD VYIVQDLMET DLYKLLKTQH LSNDHICYFL YQILRGLKYI HSANVLHRDL KPSNLLLNTT CDLKICDFGL ARVADPDHDH TGFLTEYVAT RWYRAPEIML NSKGYTKSID IWSVGCILAE MLSNRPIFPG KHYLDQLNHI LGILGSPSQE DLNCIINLKA RNYLLSLPHK NKVPWNRLFP NADSKALDLL DKMLTFNPHK RIEVEQALAH PYLEQYYDPS DEPIAEAPFK FDMELDDLPK EKLKELIFEE TARFOPGYRS

General References

Hu S., et al. (2009) Cell. 139(3):610-22. McGargill MA., et al. (2009) J Immunol. 183(8):4838-42.

DATA

SDS-PAGE



15% SDS-PAGE (3ug)

3ug by SDS-PAGE under reducing condition and visualized by coomassie blue stain.

