NKMAXBIO We support you, we believe in your research

Recombinant human PHPT1 protein

Catalog Number: ATGP0751

PRODUCT INFORMATION

Expression system

E.coli

Domain

1-125aa

UniProt No.

O9NRX4

NCBI Accession No.

NP 054891

Alternative Names

14 kDa phosphohistidine phosphatase, bA216L13.10, CGI-202, DKFZp564M173, HSPC141, PHP14, RP11-216L13.10

PRODUCT SPECIFICATION

Molecular Weight

15.9 kDa (145aa) confirmed by MALDI-TOF

Concentration

0.5mg/ml (determined by Bradford assay)

Formulation

Liquid in. 20mM Tris-HCl buffer (pH 8.0) containing 0.2M NaCl, 2mM DTT, 10% glycerol

Purity

> 95% by SDS-PAGE

Tag

His-Tag

Application

SDS-PAGE

Storage Condition

Can be stored at +2C to +8C for 1 week. For long term storage, aliquot and store at -20C to -80C. Avoid repeated freezing and thawing cycles.

BACKGROUND

Description

PHPT1 is a 125 amino acid enzyme belonging to the Janus protein family. Existing as a monomer in the cytoplasm, PHPT1 is an EDTA-insensitive phosphohistidine phosphatase. Overexpression of PHPT1 leads to specific phosphohistidine phosphatase activity towards phosphopeptide I, with no activity detected towards phosphotyrosine, phosphothreonine and phosphoserine peptides. Recombinant human PHPT1 protein, fused to His-tag at N-terminus, was expressed in E. coli and purified by using conventional chromatography techniques.



NKMAXBio We support you, we believe in your research

Recombinant human PHPT1 protein

Catalog Number: ATGP0751

Amino acid Sequence


MGSSHHHHHH SSGLVPRGSH MAVADLALIP DVDIDSDGVF KYVLIRVHSA PRSGAPAAES KEIVRGYKWA EYHADIYDKV SGDMQKQGCD CECLGGGRIS HQSQDKKIHV YGYSMAYGPA QHAISTEKIK AKYPDYEVTW ANDGY

General References

Srivastava S., et al. (2008) Proc Natl Acad Sci u S A. 105(38):14442-6. van Vlierberghe P., et al. (2006) Leukemia. 20(7):1245-53.

DATA

SDS-PAGE



15% SDS-PAGE (3ug)

3ug by SDS-PAGE under reducing condition and visualized by coomassie blue stain.

