

Recombinant human Cyclin B1 protein

Catalog Number: ATGP1612

PRODUCT INFORMATION

Expression system

E.coli

Domain

1-433aa

UniProt No.

P14635

NCBI Accession No.

NP_114172.1

Alternative Names

Cyclin B1, CCNB

PRODUCT SPECIFICATION

Molecular Weight

50.9 kDa (457aa)

Concentration

0.25mg/ml (determined by Bradford assay)

Formulation

Liquid in. 20mM Tris-HCl buffer (pH 8.0) containing 0.15M NaCl, 20% glycerol

Purity

> 85% by SDS-PAGE

Tag

His-Tag

Application

SDS-PAGE

Storage Condition

Can be stored at +2C to +8C for 1 week. For long term storage, aliquot and store at -20C to -80C. Avoid repeated freezing and thawing cycles.

BACKGROUND

Description

CCNB1, also known as cyclin B1, is a regulatory protein involved in mitosis. This protein complexes with p34 (cdc2) to form the maturation-promoting factor (MPF). Two alternative transcripts have been found, a constitutively expressed transcript and a cell cycle-regulated transcript that is expressed predominantly during G2/M phase. The different transcripts result from the use of alternate transcription initiation sites. Recombinant human CCNB1 protein, fused to His-tag at N-terminus, was expressed in E. coli and purified by using conventional chromatography techniques.

Recombinant human Cyclin B1 protein

Catalog Number: ATGP1612

Amino acid Sequence


<MGSSHHHHHH SSGLVPRGSH MGSH>MALRVT RNSKINAENK AKINMAGAKR VPTAPAATSK PGLRPRTALG
DIGNKVSEQL QAKMPMKKEA KPSATGKVID KKLKPLEKV PMLVPVPVSE PVPEPEPEPE PEPVKEEKLS PEPILVDTAS
PSPMETS GCA PAEEDLCQAF SDVILAVNDV DAEDGADPNL CSEYVKDIYA YLRQLEEEQA VRPKYLLGRE VTGNMRAILI
DWLVQVQMKF RLLQETMYMT VSIIDRFMQN NCVPKKMLQL VGV TAMFIAS KYEEMYPPEI GDFAFVTDNT YTKHQIRQME
MKILRALNFG LGRPLPLHFL RRASKIGEVD VEQHTLAKYL MELTMLDYDM VHFPPSQIAA GAFCLALKIL DNGEWTPTLQ
HYLSYTEESL LPVMQHLAKN VVMVNQGLTK HMTVKNKYAT SKHAKISTLP QLNSALVQDL AKAVAKV

General References

Antinore M J., et al. (1999). Oncogene 18 (18):2892-900

DATA

SDS-PAGE



3 μ g by SDS-PAGE under reducing condition and visualized by coomassie blue stain.