

Recombinant human TDG protein

Catalog Number: ATGP1705

PRODUCT INFORMATION

Expression system

E.coli

Domain

1-410aa

UniProt No.

Q13569

NCBI Accession No.

AAH37557

Alternative Names

Thymine-DNA glycosylase, G/T mismatch-specific thymine DNA glycosylase,

PRODUCT SPECIFICATION

Molecular Weight

48.4 kDa (433aa)

Concentration

0.5mg/ml (determined by Bradford assay)

Formulation

Liquid in. 20mM Tris-HCl buffer (pH 8.0) containing 0.4 urea, 10% glycerol, 0.1M NaCl

Purity

> 85% by SDS-PAGE

Tag

His-Tag

Application

SDS-PAGE, Denatured

Storage Condition

Can be stored at +2C to +8C for 1 week. For long term storage, aliquot and store at -20C to -80C. Avoid repeated freezing and thawing cycles.

BACKGROUND

Description

Thymine-DNA glycosylase, also known as TDG, belongs to the TDG/mug DNA glycosylase family. TDG is a nuclear protein which corrects G/T mismatches to G/C pairs by hydrolyzing the carbon-nitrogen bond between the sugar-phosphate backbone of the DNA and the mispaired thymine. TDG can also remove uracil and 5-bromouracil from mispairings with guanine. This enzyme plays a central role in cellular defense against genetic mutation caused by the spontaneous deamination of 5-methylcytosine and cytosine. Recombinant human TDG protein, fused to His-tag at N-terminus, was expressed in E. coli

Recombinant human TDG protein

Catalog Number: ATGP1705

Amino acid Sequence


<MGSSHHHHHH SSGLVPRGSH MGS>MEAENAG SYSLQQAQAF YTFPFQQLMA EAPNMAVVNE QQMPEEVPAP
APAQEPVQEA PKGRKRKPRTPTEPKQPVEPK KPVESKKS GK SAKSKEKQEK ITDTFKVKRK VDRFNGVSEA ELLTKTLPDI
LTFNLDIVII GINPGLMAAY KGHYHPGPGN HFWKCLFMSG LSEVQLNHMD DHTLPGKYGI GFTNMVERTT PGSKDLSSKE
FREGGRILVQ KLQKYQPRIA VFNGKCIYEI FSKEVFGVKV KNLEFGLQPH KIPDTETLCY GMPSSSARCA QFPRAQDKVH
YYIKLKDRLD QLKGIERNMD VQEVQYTFDL QLAQEDAKKM AVKEEKYDPG YEAAAYGGAYG ENPCSSEPCG FSSNGLIESV
ELRGESAFSG IPNGQWMTQS FTDQIPFSFN HCGTQEQQEEE SHA

General References

um S., et al. (1998) J Biol Chem. 273:20728-20736.
Hardeland u., et al. (2002) EMBO J. 21:1456-1464.

DATA

SDS-PAGE



3ug by SDS-PAGE under reducing condition and visualized by coomassie blue stain.

15% SDS-PAGE (3ug)