

Recombinant human ACAT1 protein

Catalog Number: ATGP1800

PRODUCT INFORMATION

Expression system

E.coli

Domain

34-427aa

UniProt No.

P24752

NCBI Accession No.

NP_000010

Alternative Names

acetyl-CoA acetyltransferase 1, ACAT, MAT, T2, THIL

PRODUCT SPECIFICATION

Molecular Weight

43.8 kDa (417aa) confirmed by MALDI-TOF

Concentration

1mg/ml (determined by Bradford assay)

Formulation

Liquid in. 20mM Tris-HCl buffer (pH 7.5) containing 0.1M NaCl, 10% glycerol, 1mM DTT

Purity

> 95% by SDS-PAGE

Tag

His-Tag

Application

SDS-PAGE

Storage Condition

Can be stored at +2C to +8C for 1 week. For long term storage, aliquot and store at -20C to -80C. Avoid repeated freezing and thawing cycles.

BACKGROUND

Description

ACAT1, also known as acetoacetyl-CoA thiolase, is an enzyme belongs to belongs to the membrane-bound acyltransferase family and Sterol o-acyltransferase subfamily. ACAT1 is an acetyl-CoA C-acetyltransferase enzyme. This enzyme catalyzes the reversible formation of acetoacetyl-CoA from two molecules of acetyl-CoA. It plays a role in lipoprotein assembly and dietary cholesterol absorption. In addition to its acyltransferase activity, it may act as a ligase. Recombinant human ACAT1 protein, fused to His-tag at N-terminus, was expressed in E. coli and purified by using conventional chromatography techniques.

Recombinant human ACAT1 protein

Catalog Number: ATGP1800

Amino acid Sequence


<MGSSHHHHHH SSGLVPRGSH MGS>VSKPTLK EVVIVSATRT PIGSFLGSL SLLPATKLSI AIQGAIEKAG IPKEEVKEAY
MGNVLQGGEG QAPTRQAVLG AGLPISTPCT TINKVCASGM KAIMMASQSL MCGHQDVMVA GGMESMSNVP
YVMNRGSTPY GGVKLEDLIV KDGLTDVYNK IHMGSCAENT AKKLNARNE QDAYAINSYT RSKAAWEAGK FGNEVIPVTV
TVKGGQPDVVV KEDEEYKRVD FSKVPKLKTV FQKENGTVTA ANASTLNDGA AALVLMTADA AKRLNVTPLA RIVAFADAAV
EPIDFPIAPV YAASMLKDV GLKKEDIAMW EVNEAFSLVV LANIKMLEID PQKVNINGGA VSLGHPIGMS GARIVGHLTH
ALKQGEYGLA SICNGGGGAS AMLIQKL

General References

Haapalainen, A.M., et al. (2007) *Biochemistry* 46: 4305-4321.
Thompson, S.L. et al. (1990) *J. Biol. Chem.* 265:5731-5735.

DATA

SDS-PAGE



3ug by SDS-PAGE under reducing condition and visualized by coomassie blue stain.

15% SDS-PAGE (3ug)