

Recombinant mouse MYL9 protein

Catalog Number: ATGP3077

PRODUCT INFORMATION

Expression system

E.coli

Domain

1-172aa

UniProt No.

Q9CQ19

NCBI Accession No.

NP_742116

Alternative Names

Myosin regulatory light polypeptide 9, A1327049, MLC20, Mylc2c, RLC-C, Myosin regulatory light chain 9, Myosin regulatory light chain 2 smooth muscle isoform, Myrl2

PRODUCT SPECIFICATION

Molecular Weight

22.4Da (196aa) confirmed by MALDI-TOF

Concentration

0.5mg/ml (determined by Bradford assay)

Formulation

Liquid in. Phosphate-Buffered Saline (pH 7.4) containing 10% glycerol, 1mM DTT

Purity

> 90% by SDS-PAGE

Tag

His-Tag

Application

SDS-PAGE

Storage Condition

Can be stored at +2C to +8C for 1 week. For long term storage, aliquot and store at -20C to -80C. Avoid repeated freezing and thawing cycles.

BACKGROUND

Description

MyI9 also known as myosin regulatory light polypeptide 9 is one of the numerous regulatory myosin light chains. Regulatory myosin light chains regulate contraction in smooth muscle and non-muscle cells via phosphorylation by myosin light chain kinase (MLCK). Phosphorylation of regulatory myosin light chains is catalyzed by MLCK in the presence of calcium and calmodulin and it increases the actin-activated myosin ATPase activity, thereby regulates the contractile activity. Recombinant mouse MyI9 was expressed in E. coli and purified by using

Recombinant mouse MYL9 protein

Catalog Number: ATGP3077

conventional chromatography techniques.

Amino acid Sequence


MGSSHHHHHH SSGLVPRGSH MGSMMSSKRA KAKTTKKRPQ RATSNVFAMF DQSIIQEFKE AFNMIDQNRD GFIDKEDLHD
MLASLGKNPT DEYLEGMMNE APGPINFTMF LTMFGEKLNQ TDPEDVIRNA FACFDEEASG FIHEDHLREL LTTMGDRFTD
EEVDEMYREA PIDKKGNFNY VEFTRILKHG AKDKDD

General References

Ma Z., et al. (2006) Mol. Cell. Biol. 26:9016-9034

DATA

SDS-PAGE



3ug by SDS-PAGE under reducing condition and visualized by coomassie blue stain.

15% SDS-PAGE (3ug)