

Recombinant human IL-21R protein

Catalog Number: ATGP3182

PRODUCT INFORMATION

Expression system

Baculovirus

Domain

20-232aa

UniProt No.

Q9HBE5

NCBI Accession No.

NP_851564

Alternative Names

Interleukin 21 receptor, IL-21 receptor, Novel interleukin receptor, NILR, CD360

PRODUCT SPECIFICATION

Molecular Weight

25.6 kDa (221aa)

Concentration

1mg/ml (determined by absorbance at 280nm)

Formulation

Liquid in. Phosphate-Buffered Saline (pH 7.4) containing 10% glycerol

Purity

> 90% by SDS-PAGE

Endotoxin level

< 1 EU per 1ug of protein (determined by LAL method)

Tag

His-Tag

Application

SDS-PAGE

Storage Condition

Can be stored at +2C to +8C for 1 week. For long term storage, aliquot and store at -20C to -80C. Avoid repeated freezing and thawing cycles.

BACKGROUND

Description

IL21R, also known as Interleukin-21 receptor, belongs to the type I cytokine receptors, and has been shown to form a heterodimeric receptor complex with the common gamma-chain, a receptor subunit also shared by the receptors for interleukin 2, 4, 7, 9, and 15. This receptor transduces the growth promoting signal of IL21, and is important for the proliferation and differentiation of T cells, B cells, and natural killer (NK) cells. The ligand

Recombinant human IL-21R protein

Catalog Number: ATGP3182

binding of this receptor leads to the activation of multiple downstream signaling molecules, including JAK1, JAK3, STAT1, and STAT3. Recombinant human IL21R, fused to His-tag at C-terminus, was expressed in insect cell and purified by using conventional chromatography techniques.

Amino acid Sequence

CPDLVCYTDY LQTVICILEM WNLHPSTLTL TWQDQYEELK DEATSCSLHR SAHNATHATY TCHMDVHFHM ADDIFSVNIT
DQSGNYSQEC GSFLLAESIK PAPPFNVTVT FSGQYNISWR SDYEDPAFYM LKGKLQYELQ YRNRGDPWAV SPRRKLISVD
SRSVSLLE FRKDSSYELQ VRAGPMPGSS YQGTWSEWSD PVIFQTQSEE LKEVEHHHHH H

General References


Parrish-Novak J. et al., (2000) Nature 408:57-63.

Hamming OJ. et al., (2012) J Biol Chem. 287(12):9454-60.

DATA

SDS-PAGE

(kDa)



3ug by SDS-PAGE under reducing condition and visualized by coomassie blue stain.

15% SDS-PAGE (3ug)