# NKMAXBIO We support you, we believe in your research

# Recombinant mouse IL-17A protein

Catalog Number: ATGP3280

### **PRODUCT INFORMATION**

### **Expression system**

Baculovirus

#### **Domain**

26-158aa

#### UniProt No.

062386

#### **NCBI Accession No.**

NP 034682.1

#### **Alternative Names**

Interleukin 17A, Cytotoxic T-lymphocyte-associated antigen 8, CTLA-8, Ctla-8, II17, IL-17A

# **PRODUCT SPECIFICATION**

# **Molecular Weight**

16kDa (141aa)

#### Concentration

0.5mg/ml (determined by absorbance at 280nm)

#### **Formulation**

Liquid in. Phosphate-Buffered Saline (pH 7.4) containing 10% glycerol

#### **Purity**

> 95% by SDS-PAGE

#### **Endotoxin level**

< 1 EU per 1ug of protein (determined by LAL method)

#### ıag

His-Tag

# **Application**

SDS-PAGE

# **Storage Condition**

Can be stored at +2C to +8C for 1 week. For long term storage, aliquot and store at -20C to -80C. Avoid repeated freezing and thawing cycles.

#### **BACKGROUND**

#### **Description**

IL17A, also known as Interleukin-17A, plays an important role in anti-microbial and chronic inflammation. Mature human IL17A shares 60% amino acid sequence identity with mouse and rat IL17A. It promotes protective mucosal and epidermal inflammation in response to microbial infection and induces chemokine production, neutrophil influx, and the production of antibacterial peptides. Recombinant mouse IL17A, fused to His-tag at C-



# NKMAXBio We support you, we believe in your research

# **Recombinant mouse IL-17A protein**

Catalog Number: ATGP3280

terminus, was expressed in insect cell and purified by using conventional chromatography techniques.

# **Amino acid Sequence**


AAIIPQSSAC PNTEAKDFLQ NVKVNLKVFN SLGAKVSSRR PSDYLNRSTS PWTLHRNEDP DRYPSVIWEA QCRHQRCVNA EGKLDHHMNS VLIQQEILVL KREPESCPFT FRVEKMLVGV GCTCVASIVR QAA<LEHHHHH H>

#### **General References**

Yao Z., et al. (1995) J. Immunol. 155:5483-5436. Kennedy J., et al. (1996) J. Interferon Cytokine Res. 16:611-617.

# **DATA**

#### **SDS-PAGE**



3ug by SDS-PAGE under reducing condition and visualized by coomassie blue stain.

