

# Recombinant human IL-5 protein

Catalog Number: ATGP3510

## PRODUCT INFORMATION

---

### Expression system

Baculovirus

### Domain

20-134aa

### UniProt No.

P05113

### NCBI Accession No.

NP\_000870

### Alternative Names

Interleukin-5, IL-5, EDF, TRF, T-cell replacing factor, B cell differentiation factor I, Eosinophil differentiation factor, Colony-stimulating factor eosinophil, BCGF-I

## PRODUCT SPECIFICATION

---

### Molecular Weight

14.2 kDa (124aa)

### Concentration

0.5mg/ml (determined by absorbance at 280nm)

### Formulation

Liquid in. Phosphate-Buffered Saline (pH 7.4) containing 10% glycerol

### Purity

> 90% by SDS-PAGE

### Endotoxin level

< 1 EU per 1ug of protein (determined by LAL method)

### Biological Activity

Measured in a cell proliferation assay using TF-1 human erythroleukemic cells. The ED50 range  $\leq 1.5$ ng/ml.

### Tag

His-Tag

### Application

SDS-PAGE, Bioactivity

### Storage Condition

Can be stored at +2C to +8C for 1 week. For long term storage, aliquot and store at -20C to -80C. Avoid repeated freezing and thawing cycles.

## BACKGROUND

---

# Recombinant human IL-5 protein

Catalog Number: ATGP3510

## Description

IL5, also known as interleukin-5, is a member of the hematopoietin receptor superfamily and is comprised of a cytokine-specific alpha chain and the common beta chain that is shared among these cytokines for signaling. It related hematopoietic cytokines that are important for allergic inflammation. It is a lineage-specific cytokine for eosinophilpoiesis and plays an important part in diseases associated with increased eosinophils, such as asthma. Recombinant human IL5, fused to His-tag at C-terminus, was expressed in insect cell and purified by using conventional chromatography techniques.

## Amino acid Sequence


<ADP>IPTEIPT SALVKETLAL LSTHRTLLIA NETLRIPVPV HKNHQLCTEE IFQGIGTLES QTVQGGTVER LFKNLSLIKK YIDGQKKKCG EERRRVNQFL DYLQEFLGVM NTEWIIES<HH HHHH>

## General References

Martinez-Moczygamba M., et al. (2003) J Allergy Clin Immunol. 112:653-665.  
 Milburn MV., et al. (1993) Nature. 363:172-176.


## DATA

### SDS-PAGE



3ug by SDS-PAGE under reducing condition and visualized by coomassie blue stain.

### Biological Activity



Measured in a cell proliferation assay using TF-1 human erythroleukemic cells. The ED50 range  $\leq 1.5$  ng/ml.