

Recombinant human FGF-21 (L174P) protein

Catalog Number: ATGP3789

PRODUCT INFORMATION

Expression system

E.coli

Domain

29-209aa

UniProt No.

Q9NSA1

NCBI Accession No.

NP_061986

Alternative Names

Fibroblast growth factor 21

PRODUCT SPECIFICATION

Molecular Weight

21.6 kDa (202aa) confirmed by MALDI-TOF (Molecular weight on SDS-PAGE will appear higher)

Concentration

1mg/ml (determined by Bradford assay)

Formulation

Liquid in. 20mM Tris-HCl buffer (pH 8.0) containing 10% glycerol

Purity

> 90% by SDS-PAGE

Endotoxin level

< 1 EU per 1ug of protein (determined by LAL method)

Tag

His-Tag

Application

SDS-PAGE

Storage Condition

Can be stored at +2C to +8C for 1 week. For long term storage, aliquot and store at -20C to -80C. Avoid repeated freezing and thawing cycles.

BACKGROUND

Description

FGF-21 is a unique member of FGF super family that is reported to be primarily produced by the liver. FGF-21 has recently been shown to have beneficial effects on a variety of metabolic parameters, including glucose and lipid metabolism and insulin sensitivity in animal models and is thought to exert potent antidiabetic and lipid-lowering agent of obesity and type 2 diabetes mellitus. FGF-21 activity depends on membrane protein beta

Recombinant human FGF-21 (L174P) protein

Catalog Number: ATGP3789

Klotho that physically complexes with various FGF receptors, thus conferring them the ability to bind FGF-21 and activate downstream signaling pathways. Recombinant human FGF-21 protein, fused to His-tag at N-terminus, was expressed in *E. coli* and purified by using conventional chromatography techniques.

Amino acid Sequence


<MGSSHHHHHH SSGLVPRGSH> MHPIPDSSPL LQFGGQVRQR YLYTDDAQQT EAHLEIREDG TVGGAADQSP
ESLLQLKALK PGVIQILGVK TSRFLCQRPD GALYGLSLHFD PEACSFRELL LEDGYNVYQS EAHGLPLHLP GNKSPHRDPA
PRGPARFLPL PGLPPAPPEP PGILAPQPPD VGSSDPLSMV GPSQGRSPSY AS

General References

Dostalova I., et al. (2009) *Physiol Res.* 58(1):1-7
Kharitononkov A., et al. (2008) *BioDrugs.* 22(1):37-44.

DATA

SDS-PAGE



3ug by SDS-PAGE under reducing condition and visualized by coomassie blue stain.