# NKMAXBio We support you, we believe in your research

# Recombinant mouse IL-15 protein

Catalog Number: ATGP3883

#### **PRODUCT INFORMATION**

#### **Expression system**

E.coli

#### **Domain**

49-162aa

#### UniProt No.

P48346

#### **NCBI Accession No.**

NP 032383

#### **Alternative Names**

Interleukin-15

#### **PRODUCT SPECIFICATION**

## **Molecular Weight**

17.6 kDa (152aa) confirmed by MALDI-TOF

#### Concentration

0.5mg/ml (determined by Bradford assay)

#### **Formulation**

Liquid in. Phosphate-Buffered Saline (pH 7.4) containing 10% glycerol

#### **Purity**

> 90% by SDS-PAGE

#### **Endotoxin level**

< 1 EU per 1ug of protein (determined by LAL method)

#### **Biological Activity**

Measured in a cell proliferation assay using CTLL2 mouse cytotoxic T cells. The ED50 for this effect is less or equal to 20ng/ml.

## **Tag**

His-Tag

### **Application**

SDS-PAGE, Bioactivity

#### **Storage Condition**

Can be stored at +2C to +8C for 1 week. For long term storage, aliquot and store at -20C to -80C. Avoid repeated freezing and thawing cycles.

### **BACKGROUND**



# NKMAXBio we support you, we believe in your research

# Recombinant mouse IL-15 protein

Catalog Number: ATGP3883

## **Description**

Interleukin (IL) -15 is a pleiotropic cytokine that plays a pivotal role in both innate and adaptive immunity. IL-15 is unique among cytokines due to its participation in a trans signaling mechanism in which IL-15 receptor alpha (IL-15R alpha) from one subset of cells presents IL-15 to neighboring IL-2R beta/gamma (c) -expressing cells. This newly discovered cytokine is produced by a wide variety of cells and tissues and like IL-2, stimulates the proliferation of T-lymphocytes and induces the generation of cytotoxicT-lymphocytes (CTLs) and lymphokine-activated killer (LAK) cells. Recombinant Mouse IL-15 was expressed in E. coli and purified by using conventional chromatography techniques.

## **Amino acid Sequence**


<MRGSHHHHHH GMASMTGGQQ MGRDLYDDDD KDRWGSHM>NW IDVRYDLEKI ESLIQSIHID TTLYTDSDFHPSCKVTAMNC FLLELQVILH EYSNMTLNET VRNVLYLANS TLSSNKNVAE SGCKECEELE EKTFTEFLQS FIRIVQMFIN TS

#### **General References**

J G Giri., et al. (1994) EMBO. 13(12): 2822-2830. Lorenzen I., et al. (2006). J Biol Chem. 281(10):6642-7.

### **DATA**

#### **SDS-PAGE**




3ug by SDS-PAGE under reducing condition and visualized by coomassie blue stain.

# **Biological Activity**

# NKMAXBio We support you, we believe in your research

# **Recombinant mouse IL-15 protein**

Catalog Number: ATGP3883



Mouse IL-15 stimulates cell proliferation of the CTLL2 mouse cytotoxic T cells. The ED50 range is less or equal to 20 ng/ml.

