

TEST003

Catalog Number: ATGA11111

PRODUCT INFORMATION

Catalog number

ATGA11111

Clone No.

dsddd

Product type

ddd

UnitProt No.

UniProt No.

NCBI Accession No.

1111

Alternative Names

111

PRODUCT SPECIFICATION

Antibody Host

Rabbit

Reacts With

Human

Concentration

1ml/mg

Formulation

dddd

Immunogen

dddd

Isotype

IgG1 kappa

Purification Note

By protein-A affinity chromatography

Application

WB, ICC/IF, IHC, FACS

Usage

vvv

Storage

2-8

TEST003

Catalog Number: ATGA11111

BACKGROUND

Description


111

General References

dddd


DATA

Western blot analysis (WB)



WB Comment and Pic

Immunocytochemistry/Immunofluorescence (ICC/IF)



IF Comment And Pic

Flow cytometry (FACS)

TEST003

Catalog Number: ATGA11111

FACS Comment and Pic

