

PRODUCT INFORMATION

Catalog number

TEST001

Clone No.

44

Product type

UnitProt No.

NCBI Accession No.

111

Alternative Names

111

PRODUCT SPECIFICATION

Antibody Host

Mouse

Reacts With

Human

Concentration

11

Formulation

11

Immunogen

Isotype

IgG1 kappa

Purification Note

By protein-A affinity chromatography

Application

SDS-PAGE

Usage

44

Storage

11

BACKGROUND

Description




111


General References

11

DATA


Western blot analysis (WB)






520x520


Immunocytochemistry/Immunofluorescence (ICC/IF)



400x400




520x520



400x400


400x400 size/520x520

Catalog Number: TEST001




520x520

Flow cytometry (FACS)




400x400



520x520

Additional Information



11