

Recombinant human DC-SIGN/CD209 protein

Catalog Number: ATGP4010

PRODUCT INFORMATION

Expression system

Baculovirus

Domain

60-404aa

UniProt No.

Q9NNX6

NCBI Accession No.

NP_066978

Alternative Names

CDSIGN, CLEC4L, CD209 antigen isoform 1, CD209 antigen, C-type lectin domain family 4 member L, Dendritic cell-specific ICAM-3-grabbing non-integrin 1, DC-SIGN, DC-SIGN1, CD_antigen, CD209

PRODUCT SPECIFICATION

Molecular Weight

66.5kDa (587aa)

Concentration

0.25mg/ml (determined by Absorbance at 280nm)

Formulation

Liquid. In Phosphate-Buffered Saline (pH 7.4) containing 10% glycerol

Purity

> 90% by SDS - PAGE

Endotoxin level

< 1 EU per 1ug of protein (determined by LAL method)

Tag

hIgG-His-Tag

Application

SDS-PAGE

Storage Condition

Can be stored at +2C to +8C for 1 week. For long term storage, aliquot and store at -20C to -80C. Avoid repeated freezing and thawing cycles.

BACKGROUND

Description

DC-SIGN/CD209, also known as Dendritic Cell-specific ICAM-3 Grabbing Non-integrin, is a member of the C-type lectin family. It is Pathogen-recognition receptor expressed on the surface of immature dendritic cells (DCs) and involved in initiation of primary immune response. DC-SIGN on macrophages recognises and binds with high

Recombinant human DC-SIGN/CD209 protein

Catalog Number: ATGP4010

affinity to high-mannose type N-glycans, a class of pathogen associated molecular patterns PAMPs commonly found on viruses, bacteria and fungi. This binding interaction activates phagocytosis. It has a high affinity for the ICAM3 molecule and binds various microorganisms by recognizing high-mannose-containing glycoproteins on their envelopes and especially functions as receptor for several viruses such as HIV and Hepatitis C. Besides functioning as an adhesion molecule, recent study has also shown that DC-SIGN can initiate innate immunity by modulating toll-like receptors, though the detailed mechanism is not yet known. It together with other C-type lectins is involved in recognition of tumors by dendritic cells. It is also a potential engineering target for dendritic cell based cancer vaccine. Recombinant human DC-SIGN/CD209, fused to hIgG-His-tag at C-terminus, was expressed in insect cell and purified by using conventional chromatography techniques.

Amino acid Sequence


<ADP>VSKVPSS ISQEQSRQDA IYQNLTLKA AVGELSEKSK LQEIQELTQ LKAAVGELPE KSKLQEIYQE LTRLKAAVGE LPEKSKLQEI YQELTWLAAA VGELPEKSKM QEIQELTRL KAAVGELPEK SKQQEIYQEL TRLKAAVGE PEKSKQQEIY QELTRLKAAV GELPEKSKQQ EIYQELTLK AAVERLCHPC PWEWTFQGN CYFMSNSQRN WHDSITACKE VQAQLVVIKS AEEQNFLQLQ SSRNRFWTM GLSDLNQEGT WQWVDGSPLL PSFKQYWNRG EPNNVGEEDC AEFSGNGWND DKCNLAKFWI CKKSAASCSR DEEQFLSPAP ATPNPPPA<VE PKSCDKTHTC PPCAPELLG GPSVFLFPPK PKDTLMISRT PEVTCVVVDV SHEDPEVKFN WYVDGVEVHN AKTKPREEQY NSTYRVVSVL TVLHQDWLNG KEYKCKVSNK ALPAPIEKTI SKAKGQPREP QVYTLPPSRD ELTKNQVSLT CLVKGFYPSD IAVEWESNGQ PENNYKTPP VLDSGDGFFL YSKLTVDKSR WQQGNVFSCS VMHEALHNHY TQKSLSLSPG KHHHHHH>

General References

Liu, W. et al, (2004) J. Biol. Chem. 279:18748.
 Anthony, R.M. et al, (2008) Proc. Natl. Acad. Sci. USA 105:19571.

DATA

SDS-PAGE



3ug by SDS-PAGE under reducing condition and visualized by coomassie blue stain