NKMAXBIO We support you, we believe in your research

Recombinant human MKK6 protein

Catalog Number: ATGP4041

PRODUCT INFORMATION

Expression system

Baculovirus

Domain

1-334aa

UniProt No.

P52564

NCBI Accession No.

NP 002749

Alternative Names

Dual specificity mitogen-activated protein kinase kinase 6, MAPKK6, MEK6, MKK6, PRKMK6, SAPKK-3, SAPKK3, Dual specificity mitogen-activated protein kinase kinase 6 isoform1, MAP2K6, MAP kinase kinase 6, MAPK/ERK kinase 6, SAPK kinase 3, SKK3

PRODUCT SPECIFICATION

Molecular Weight

38.3 kDa (340aa) confirmed by MALDI-TOF

Concentration

1mg/ml (determined by absorbance at 280nm)

Formulation

Liquid. In Phosphate-Buffered Saline (pH 7.4) containing 20% glycerol

Purity

> 90% by SDS-PAGE

Endotoxin level

< 1 EU per 1ug of protein (determined by LAL method)

Tag

His-Tag

Application

SDS-PAGE

Storage Condition

Can be stored at +2C to +8C for 1 week. For long term storage, aliquot and store at -20C to -80C. Avoid repeated freezing and thawing cycles.

BACKGROUND

Description

MKK6, also known as dual specificity mitogen-activated protein kinase kinase 6, is a member of the dual specificity protein kinase family, which functions as a mitogen-activated protein (MAP) kinase kinase. MAP



NKMAXBIO We support you, we believe in your research

Recombinant human MKK6 protein

Catalog Number: ATGP4041

kinases, also known as extracellular signal-regulated kinases (ERKs), act as an integration point for multiple biochemical signals. This protein phosphorylates and activates p38 MAP kinase in response to inflammatory cytokines or environmental stress. As an essential component of p38 MAP kinase mediated signal transduction pathway, this gene is involved in many cellular processes such as stress-induced cell cycle arrest, transcription activation and apoptosis. Recombinant human MKK6, fused to His-tag at C-terminus, was expressed in insect cell and purified by using conventional chromatography techniques.

Amino acid Sequence


MSQSKGKKRN PGLKIPKEAF EQPQTSSTPP RDLDSKACIS IGNQNFEVKA DDLEPIMELG RGAYGVVEKM RHVPSGQIMA VKRIRATVNS QEQKRLLMDL DISMRTVDCP FTVTFYGALF REGDVWICME LMDTSLDKFY KQVIDKGQTI PEDILGKIAV SIVKALEHLH SKLSVIHRDV KPSNVLINAL GQVKMCDFGI SGYLVDSVAK TIDAGCKPYM APERINPELN QKGYSVKSDI WSLGITMIEL AILRFPYDSW GTPFQQLKQV VEEPSPQLPA DKFSAEFVDF TSQCLKKNSK ERPTYPELMQ HPFFTLHESK GTDVASFVKL ILGD

General References

Chen Z, et al (2001) J Biol Chem. 276(19):16070-16075. Raingeaud J, et al (1996) Mol Cell Biol. 16:1247-1255.

DATA

SDS-PAGE



3ug by SDS-PAGE under reducing condition and visualized by coomassie blue stain

