NKMAXBIO We support you, we believe in your research

Recombinant human IL-17RB protein

Catalog Number: ATGP4039

PRODUCT INFORMATION

Expression system

Baculovirus

Domain

18-292aa

UniProt No.

O9NRM6

NCBI Accession No.

NP 061195.2

Alternative Names

Interleukin-17 receptor B, CRL4, EVI27, IL17BR, IL17RH1, IL-17 receptor B, IL-17RB, Cytokine receptor-like 4, IL-17 receptor homolog 1, Interleukin-17B receptor

PRODUCT SPECIFICATION

Molecular Weight

57.3 kDa (514aa) Reducing conditions

Concentration

0.25mg/ml (determined by Bradford assay)

Formulation

Liquid in. Phosphate-Buffered Saline (pH 7.4) containing 10% glycerol

Purity

> 85% by SDS-PAGE

Endotoxin level

< 1 EU per 1ug of protein (determined by LAL method)

Biological Activity

Measured by its binding ability in a functional ELISA with Human IL25 (CAT# ATGP4018). The ED50 range \leq 2 ug/ml.

Tag

hlgG-His-Tag

Application

SDS-PAGE, Bioactivity

Storage Condition

Can be stored at +2C to +8C for 1 week. For long term storage, aliquot and store at -20C to -80C. Avoid repeated freezing and thawing cycles.



NKMAXBIO We support you, we believe in your research

Recombinant human IL-17RB protein

Catalog Number: ATGP4039

BACKGROUND

Description

IL-17RB, also known as Interleukin-17 receptor B, is a cytokine receptor. This receptor specifically binds to IL17B and IL17E (IL25), but does not bind to IL17A or IL17C. It has been shown to mediate the activation of NFkappaB and the production of IL8 induced by IL17E. It was found to be significantly up-regulated during intestinal inflammation, which suggested the immunoregulatory activity of this receptor. Recombinant human IL-17RB, fused to hIgG-His-tag at C-terminus, was expressed in insect cell and purified by using conventional chromatography techniques.

Amino acid Sequence


REPTVQCGSE TGPSPEWMLQ HDLIPGDLRD LRVEPVTTSV ATGDYSILMN VSWVLRADAS IRLLKATKIC VTGKSNFQSY SCVRCNYTEA FQTQTRPSGG KWTFSYIGFP VELNTVYFIG AHNIPNANMN EDGPSMSVNF TSPGCLDHIM KYKKKCVKAG SLWDPNITAC KKNEETVEVN FTTTPLGNRY MALIQHSTII GFSQVFEPHQ KKQTRASVVI PVTGDSEGAT VQLTPYFPTC GSDCIRHKGT VVLCPQTGVP FPLDNNKSKP GGWLP<LEPKS CDKTHTCPPC PAPELLGGPS VFLFPPKPKD TLMISRTPEV TCVVVDVSHE DPEVKFNWYV DGVEVHNAKT KPREEQYNST YRVVSVLTVL HQDWLNGKEY KCKVSNKALP APIEKTISKA KGQPREPQVY TLPPSRDELT KNQVSLTCLV KGFYPSDIAV EWESNGQPEN NYKTTPPVLD SDGSFFLYSK LTVDKSRWQQ GNVFSCSVMH EALHNHYTQK SLSLSPGKHH HHHH>

General References

Ma XJ,. et al. (2006) J Clin Oncol. 24:4611-4619. Schwandner R,. et al. (2000) J. Exp. Med. 191:1233-1240.

DATA

SDS-PAGE




3ug by SDS-PAGE under reducing condition and visualized by coomassie blue stain

Biological Activity



Recombinant human IL-17RB protein

Catalog Number: ATGP4039



Human IL25 (CAT# ATGP4018) is coated at 2 ug/ml (100 ul/well) can bind Human IL17RB. The ED50 range \leq 2 ug/ml.

