

Mouse Progesterone receptor ELISA Kit

Catalog Number: ATGE0101

PRODUCT INFORMATION

Catalog number


ATGE0101

PRODUCT SPECIFICATION

Description

Cat.No E0490Mo Standard Curve Range: 0.5-150ng/ml Sensitivity: 0.26ng/ml Size: 96 wells / 48 wells Storage: Store the reagents at 2-8°C. For over 6-month storage refer to the expiration date keep it at -20°C. Avoid repeated thaw cycles. If individual reagents are opened it is recommended that the kit be used within 1 month. *This product is for research use only, not for use in diagnosis procedures. It's highly recommended to read this instruction entirely before use. Precision Intra-Assay Precision (Precision within an assay) Three samples of known concentration were tested on one plate to assess intra-assay precision. Inter-Assay Precision (Precision between assays) Three samples of known concentration were tested in separate assays to assess inter-assay precision. $CV(\%) = SD/mean \times 100$ Intra-Assay: CV Intended Use This sandwich kit is for the accurate quantitative detection of Mouse Progesterone receptor (also known as PGR) in serum, plasma, cell culture supernates, Ascites, tissue homogenates or other biological fluids.

SDS-PAGE



Storage

2-8C

Product type